


Jan 13, 2025

3. Data Generation - Systems with an Area Detector

 In 2 collections

DOI

dx.doi.org/10.17504/protocols.io.bp2l62eorgqe/v1

Mathias Hammer¹, Britta Schroth-Diez², Luis-Francisco Acevedo-Hueso³, Gert-Jan Bakker⁴, Laszlo Barna⁵, Sebastian Beer⁶, Valeria Berno⁷, Yann Cesbron⁸, Orestis Faklaris⁹, Nathalie Gaudreault¹⁰, Thomas Guilbert¹¹, Rainer Heintzmann¹², Gabriel G Martins¹³, David McFadden¹⁴, Glyn Nelson¹⁵, Roland Nitschke¹⁶, Kees van der Oord¹⁷, Sathya Srinivasan¹⁸, Andre Zeug¹⁹, David Grunwald²⁰

¹UMASS Chan Medical School, RNA Therapeutics Institute, Worcester, MA, USA;

²Max Planck Institute of Molecular Cell Biology and Genetics, Dresden, Germany;

³Neuro-AI Monitoring and Quantum Medical Diagnostics;

⁴Radboud University Medical Center, Medical BioSciences department, Nijmegen, the Netherlands;

⁵Institute of Experimental Medicine, Budapest; ⁶Hamamatsu Photonics Deutschland GmbH;

⁷ALEMBIC Advanced Light and Electron Microscopy BiImaging Center, San Raffaele Scientific Institute;

⁸Institute of Science and Technology Austria;

⁹MRI, Biocampus, University of Montpellier, CNRS, INSERM, Montpellier, France;

¹⁰Allen Institute for Cell Science, Seattle, WA, USA;

¹¹Institut Cochin, INSERM U1016, CNRS UMR 8104, Université Paris Cité, Paris, France;

¹²Leibniz Institute of Photonic Technology; ¹³Gulbenkian Institute for Molecular Medicine;

¹⁴Friedrich Schiller University Jena; ¹⁵Newcastle University, BiImaging Unit, UK;

¹⁶Life Imaging Center, University of Freiburg, Germany; ¹⁷Nikon Europe BV; ¹⁸ONPRC;

¹⁹Cellular Neurophysiology, Hannover Medical School, Hannover, Germany;

²⁰UMass Chan Medical School, RNA Therapeutics Institute, Worcester, MA, USA

QUAREP-LiMi

Tech. support email: info@quarep.org

Britta Schroth-Diez

Max Planck Institute of Molecular Cell Biology and Genetics,...

OPEN  ACCESSDOI: dx.doi.org/10.17504/protocols.io.bp2l62eorgqe/v1External link: <https://quarep.org/>

Protocol Citation: Mathias Hammer, Britta Schroth-Diez, Luis-Francisco Acevedo-Hueso, Gert-Jan Bakker, Laszlo Barna, Sebastian Beer, Valeria Berno, Yann Cesbron, Orestis Faklaris, Nathalie Gaudreault, Thomas Guilbert, Rainer Heintzmann, Gabriel G Martins, David McFadden, Glyn Nelson, Roland Nitschke, Kees van der Oord, Sathya Srinivasan, Andre Zeug, David Grunwald 2025. 3. Data Generation - Systems with an Area Detector. [protocols.io](https://www.protocols.io) <https://dx.doi.org/10.17504/protocols.io.bp2l62eorgqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: September 08, 2024

Last Modified: January 13, 2025

Protocol Integer ID: 107124

Keywords: detection system, detector, camera, area detector, characterization, in-homogeneous, acquisition, measurement, CMOS, sCMOS, CCD, EMCCD, dynamic range, noise

Disclaimer

This protocol collection was developed by members of **WG2 "Detection System Performance"** of the international consortium **QUAREP-LiMi**.

The Consortium for Quality Assessment and Reproducibility for Instruments and Images in Light Microscopy (QUAREP-LiMi), formed by the global community of practitioners, researchers, developers, service providers, funders, publishers, policy makers and industry related to the use of light microscopy, is committed to democratizing access to quantitative and reproducible light microscopy and the data generated by it.

This protocol collection has undergone the internal approval process of QUAREP-LiMi.

It is not the intention of this protocol collection to supplant the need for independent professional judgement, advice, diagnosis or treatment. Any action taken or withheld on the basis of the information presented here is undertaken at the user's own risk. The user agrees that neither the company nor any of the authors, contributors, administrators, or other individuals associated with protocols.io or QUAREP-LiMi can be held liable for any injuries, damages or losses incurred as a result of the user's use of the information contained in or linked to this protocol or any of our websites, applications, or services.

Abstract

This protocol describes the measurement procedure to produce dark and inhomogeneously illuminated images with a light microscope system equipped with area detectors. The protocol can be followed according to the three different aims of microscope characterization described in the introduction of the protocol collection "Characterization of the Photon Conversion Factor, Noise, and Dynamic Range of Light Microscope Detection Systems". The protocol uses the microscope slide for inhomogeneous illumination described in protocol 2 and generates data which can be analyzed with protocol 5, to obtain the photon conversion factor, readnoise, and dynamic range of the detection system.

Guidelines

Please refer to protocol "**1. Introduction - Background and Aims**" of this collection for a more detailed description of a detection system, its parameters and various aims for performance monitoring.

This protocol is part of a collection of protocols developed by QUAREP-LiMi WG2 for characterising detection system performance.

We recommend to perform the described measurements one to two times a year, or if you notice unusual behaviour of your detector.

Materials

High dynamic range fluorescent test sample. E.g., marker pen slide as discribed in **protocol 2** of this collection.

Safety warnings

- ! Ensure you follow general lab safety guidelines for radiation sources and chemicals as outlined within your organisation.



Laser safety and regulations

- Please refer to the documentation provided by the manufacturer for additional warnings and preventive, protective equipment (PPE) requirements (e.g. laser safety goggles). Always consult your local Laser Safety Officer or Radiation Safety Officer and refer to your laboratory safety documentation for more information.
- You can also consult your Laser Safety Standards ANSI Z136 in North America, SUVA 66049.D in Europe, and BS EN 60825-1 in the UK. Additionally, laser safety standards and regulations are covered by IEC norm 60825-1, and LED eye safety standards and regulations are covered by IEC norm 62471 in Europe.

Safety information

Hazardous, visible, or invisible radiation from lasers, lamps, and other light sources used for microscopy **can cause permanent damage to the retina, skin burns, and fire**. Always follow proper laser safety protocols for your equipment and situation.

Before start

This protocol can be used to determine experiment QC (Aim 1), instrument QC (Aim 2) or characterisation of the system (Aim 3). In each case, the notes describe the breakdown for each Aim. For further explanations of each aim, please refer to protocol "**1. Introduction - Background and Aims**" of this collection.

Make sure the microscope is turned on at least one hour prior to the measurements to allow the system to stabilize. All measurements in this protocol must be performed without any ambient light contributing to the signal.

For consistency, remember to always record your settings with your images. All settings must be kept unchanged to be able to reliably monitor system performance over time.

Set up the microscope

- 1 Select the appropriate objective.

Note

For **Aim 1 - experiment QC**: use the objective according to the settings of your experiment

For **Aim 2 - instrument QC and Aim 3 - system characterization**: use a 10x objective (or 20x if a 10x is not available); if measurements are done repeatedly, always use the same objective.

- 2 Install a high-dynamic range test sample on the stage.

Note

It is essential to have both, very bright as well as very dark, areas in the sample to ensure a high dynamic range signal distribution. A maximum contrast (signal-to-background) ratio of at least 100:1 is ideal.

Suitable test samples can be prepared following our [fluorescent test sample protocol](#).

- 3 Focus and position your test sample.

Note

For the blurry edge test sample (prepared according to the protocol mentioned above), bring the test sample in focus and adjust its position such that the bright area fills about one-third of the field of view.




Image and histogram of the test sample (black: linear scaling; grey: logarithmic scaling).

- 4 Defocus the sample slightly to create a continuous intensity distribution.

Note

Defocussing by about one micrometer creates a blurry edge.




Image and histogram of the sample after defocussing. Histogram: linear (black) and logarithmic (grey) scaling. The intensity distribution in the histogram is flattened out in comparison to in focus images.

Adjust the detector settings

- 5 Adjust the acquisition parameters according to your requirements.

Note

Acquisition parameters include: detector settings such as exposure time, gain settings, temperature set point, read out rates, bit depth, binning, as well as other parameters, like the illumination power of the light source or the filters used in the optical path of the detection system.

Note

For **Aim 1**: Use your experiment settings.
For **Aim 2**: Use always the same analog gain setting, the same operational mode, and the same exposure times (between 50 and 100ms).
For **Aim 3**: We recommend to use all standard analog gain settings, standard operational modes, and an exposure time between 50 and 100ms.
For further considerations please see step 6 and respective sub steps.

Note

Check any other settings that may affect the image or image integrity. Pre-processing filters (such as pixel defect removal, camera autocontrast, or autodeblur) will affect the calibration results. We recommend turning these off for calibration and quality control to ensure repeatability and consistency, especially for Aims 2 and 3.

- 6 Choose between saturation-level calibration and arbitrary-level calibration. Adjust your acquisition settings accordingly.

Note**1. Saturation-level calibration:**

The calibration is performed at intensities ranging from very low up to the saturation level, covering the entire dynamic range of the detector. A saturation-level calibration can be useful as a repeatable "one size fits all" approach. Saturation-level calibration will, in addition, yield information on the total dynamic range of the detector.

2. Arbitrary-level calibration:

The calibration is performed at an arbitrary intensity level below the saturation level. Some experiments benefit from a calibration performed close to the intensities of the actual experiment since detectors often show significant performance differences depending on the charge load of its detection elements. Arbitrary-level calibration yields information about the total number of detected photons.


- 6.1 **For saturation-level calibration:** Adjust the parameters so that around 1% of the pixels are saturated.

Note

The precise number of saturated pixels is not important.

The percentage of saturated pixels can be determined after checking the results according to the analysis protocol.

Please be aware that oversaturation may harm the detector. This is especially true for EMCCDs.




A saturated image and its histogram (black: linear intensity distribution; grey: logarithmic intensity distribution). Highlighted in red are the saturated pixels.

Note

For EM-CCDs grayscale values corresponding to saturation greatly depend on the camera settings. Detector saturation might occur before the A/D conversion. Saturation-level calibration is a good tool for determining this.

- 6.2 **For arbitrary-level calibration:** If you want to calibrate for an illumination range that is less than 50% of the detector saturation level, adjust the acquisition parameters so that the brightest pixels in the image are 20% higher than your target maximum intensity.

Note

This figure shows a non-saturated image and its histogram (black: linear intensity distribution; grey: logarithmic intensity distribution). There are no saturated pixels. The histogram spans across about a third of the 16bit range.

Data acquisition

7 Acquire a bright image series and a dark image series.

7.1 **Bright image series:** Take a series of at least **20 images** with the selected settings and with the respective illumination turned **ON**.

Note

Some detectors show noticeably different characteristics during a "run-in" period caused by for instance changes in offset or sensor temperature. The "run-in" time corresponds to the time until the intensity between images stays constant, apart from noise.

We suggest a series of at least 20 images for the initial measurement. Adjust the number of images to the "run-in" time plus 20 images, remove the "run-in" images from the series before the analysis.

Note

Detectors can perform differently depending on the acquisition conditions (exposure time, temperature and illumination level). If you wish to calibrate multiple conditions, you will need to acquire a new series for each set of conditions (**Aims 1 through 3**).

7.2 **Dark image series:** Take a series of at least **20 images** with the same respective settings, but with the illumination turned **OFF**.

Note

Dark images are necessary for determining a detector's offset value and dark noise.

Note

To ensure that there is no ambient light contributing to the signal, either the full detection light path needs to be blocked or covered, or the entire microscopy room must be darkened.

Note

If the detector settings are kept the same and only illumination or optical path settings are changed, the same 20 dark images can be used for analysis. No additional dark image series has to be acquired in this case.

- 8 Record all acquisition parameters (metadata) with each image series.

Note

Metadata are needed for further data handling and analysis and thus must be kept together with the image data. We recommend to store data in the microscope manufacturer's data format (where available). In most cases this will ensure storage of all the required metadata.

- 9 Repeat steps [⇒ go to step #5](#) through [⇒ go to step #8](#) for each additional set of acquisition parameters or illumination levels for which you want to perform a calibration.

Analysis

- 10 For respective analysis steps, please refer to the [Analysis protocol](#).

Acknowledgements

 Funding acknowledgement.pdf 50KB

 CRediT_author_contributions_3_Dat...



protocols.io